NN31545.1863

BIBLIOTHEEK
STARINGGEBOUW

ICW note 1863
April 1988

I

nota

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METHODS FOR THE DETERMINATION OF CHEMICAL COMPOSITION
OF WATER SAMPLES FROM COLUMN EXPERIMENTS

English and Indonesian version

A. v.d. Toorn, ir. K. Nugroho, drs. J. Harmsen

0000 0335 6355

Nota's (Notes) of the Institute are a means of internal communication and not a publication. As such their contents vary strongly, from a simple presentation of data to a discussion of preliminary research results with tentative conclusions. Some notes are confidential and not available to third parties if indicated as such

1912-291699* 20 JULI 1988

C O N T E N T S

	page
1. INTRODUCTION	1
2. MAINTENANCE OF INSTRUMENTATION	3
2.1. General	3
2.2. Daily	3
2.3. Weekly	3
2.4. Monthly	3
2.5. Maintenance of electrodes	4
3. DISSOLVED OXYGEN AND PH DETERMINATION	5
3.1. OXI 91 Oxygenmeter	5
3.2. Measuring of pH (pH 535 multical)	6
4. CONDUCTIVITY	9
4.1. WTW Laboratory meter	9
4.2. Field equipment	10
4.3. Chemicals	12
5. DETERMINATION OF CHLORIDE CONTENT USING ELECTRODE	13
5.1. Water samples	13
5.2. Calculation	13
5.3. Chemicals	13
6. DETERMINATION OF BICARBONATE AND CARBONATE CONTENT	14
6.1. Water samples	14
6.2. Calculation	14
6.3. Chemicals	14
7. PHOTOMETRIC DETERMINATION ON NITRATE CONTENT	15
7.1. Apparatus	15
7.2. Sample	15
7.3. Procedure	15
7.4. Calculations	16
7.5. Chemicals	16
7.6. References	17

	page
8. PHOTOMETRIC DETERMINATION OF SULPHATE CONTENT	18
8.1. Instrumentation	18
8.2. Reagents	18
8.3. Procedure	18
8.4. Calculation	19
8.5. Chemicals/reagents	20
9. DETERMINATION OF SODIUM AND POTASSIUM WITH FLAMEPHOTOMETER	21
9.1. Calibration of the meter (monthly)	21
9.2. Measurements	21
9.3. Aspirate	21
9.4. Calculation	22
9.5. Chemicals	22
10. DETERMINATION OF NA, K, CA, MG, FE AND MN WITH THE ATOMIC ABSORPTION SPECTROPHOTOMETER	23
11. PHOTOMETRIC DETERMINATION OF ALUMINIUM CONTENT	27
11.1. Apparatus	27
11.2. Samples	27
11.3. Procedure	27
11.4. Calculation	28
11.5. Chemicals	29
12. ADJUSTING NORMALITIES	30
12.1. Normality of AgNO_3	30
12.2. Normality of HCl	30
13. CHLORIDE DETERMINATION USING INDICATOR	32
13.1. Watersamples	32
13.2. Calculation	32
13.3. Chemicals	32
14. BICARBONATE AND CARBONATE DETERMINATION USING TITRATION METHODS WITH INDICATOR	33
14.1. Watersamples	33
14.2. Calculation	33
14.3. Chemicals	33
15. CALCIUM AND MAGNESIUM DETERMINATION USING TITRATION	34
15.1. Calcium and magnesium (total hardness)	34
15.2. Calcium determination	35
15.3. Magnesium determination	35
15.4. Chemicals	36
LITERATURE	39

1. INTRODUCTION

In the chemical laboratory, water samples from column experiments have to be analysed for several parameters as presented in table 1.

Modern equipment like spectrophotometers, flamephotometer, atomic absorption spectrophotometers, digital pH and mv meters, conductivity meters will be used. Also digital burets, an analytical balance, a freeze dryer, and two apparatus to supply distilled water will be available. The instruments used for a certain determination are also given in table 1. The last column in the table gives the principle of the method used. Some parameters can be analysed with different methods, which is important when one fails.

Table 1. Methods of the chemical laboratory

Parameter	Instrument	Principle	Related chapter
pH	pH meter	electrode	3
Conductivity	conductivity meter	electrode	4
Chloride	mV meter	titration/electrode	5
Carbonate	pH meter	titration/electrode	6
Bicarbonate	pH-meter	titration/electrode	6
Nitrate	spectrophotometer	color development	7
Sulphate	spectrophotometer	color development	8
Sodium	flamephotometer	atomic emission	9
Potassium	flamephotometer	atomic emission	9
Magnesium	AAS	atomic absorbtion	10
Calcium	AAS	atomic absorbtion	10
Iron	AAS	atomic absorbtion	10
Manganese	AAS	atomic absorbtion	10
Aluminium	spectrophotometer	color development	11

Parameter	Alternative way	Principle	Related chapter
pH	pH paper	color development	-
Chloride	digital buret	titration/color	13
Carbonate	digital buret	titration/color	14
Bicarbonate	digital buret	titration/color	14
Sodium	AAS	atomic absorbtion	10
Potassium	AAS	atomic absorbtion	10
Calcium	digital buret	titration/color	15
Magnesium	digital buret	titration/color	15
Nitrate	colorimetric*	color development	-
Sulphate	colorimetric*	color development	-
Aluminium	colorimetric*	color development	-

* Visual comparison with the standard series.

This is a less accurate method.

2. MAINTENANCE OF INSTRUMENTATION

2.1. General

To work in a proper way with laboratory instrumentation, equipment needs maintenance. This can be a simple cleaning, but it can also be necessary to disassemble an instrument. Especially because of the problem of dust maintenance is very important. Without this every instrument will fail after a short period. Maintenance can be done daily to yearly.

2.2. Daily

- Place electrodes in distilled water or in another solution which is recommended.
- Clean all tubes, beakers, pipettes etc which has been used.
- Check if all gas bottles are closed.
- If there are samples, calibrate the conductivity and the pH-meter.

2.3. Weekly

- Check level of reference electrodes and if necessary refill with the proper solution.
- Renew the solution in which the electrodes are stored.
- Clean balances.
- Check air filter of the atomic absorption and remove condensed water. If necessary replace the filter.
- Suck some distilled water through the flowcell of the Spectronic 20.

2.4. Monthly

- Adjust normalities of all titration solutions and compare with previous ones.

2.5. Maintenance of electrodes

2.5.1. Silverelectrode

When AgCl has precipitated on the electrode, it slowly reacts. Remove this AgCl with a tissue.

2.5.2. Saturated K_2SO_4 reference electrode

This electrode must always be filled with saturated K_2SO_4 up to 1 cm under the filling hole. If the electrode is not in use, it has to be stored in a saturated K_2SO_4 solution. If the signal to the mV meter is not stable it is possible that the membrane is clogged with AgCl. This can be removed with a piece of fine abrasive sandpaper.

3. DISSOLVED OXYGEN AND PH DETERMINATION

3.1. OXI 91 Oxygen-meter

- Plug in Oxy-genelectrode at the back side of the instrument (figure 1).


Fig. 1. OXI 91 oxygen-meter

- Main switch at ' O_2 ' and the upper switch at 'MBAR'.
- Read airpressure at the airpressure meter and set the display of the oxygen-meter at the right value, with the left knob at the front side of the instrument.
- Upper switch to '% SAT'.
- Pull the black stick out of the Oxical calibration cell and put the oxygen-electrode in (take care that the electrode is dry).
- Wait a few minutes.

- Set the oxygen-meter with the right knob at the front to 102%.
- Take the electrode out of the calibration cell and put the black stick in the cell again.
- The instrument is ready now for measurements.
- Measure possibilities '% SAT' and mg/l.

3.2. Measuring of pH (pH 535 multical)

3.2.1. General

- Stir all solutions with a magnetic stirrer.
- Rinse the electrodes with distilled water before use in another solution.

If not in use keep the electrodes in a 3 m KCl solution.

3.2.2. Calibration of the pH meter

- Turn on the instrument with 'MAINS' at backside of the instrument.


Fig. 2. Front figure of pH meter 535 multical


Fig. 3. Back side of pH meter 535 multical

- Press button 'MODE/RUN' till function 'CAL PH' appears on the display.
 - The pH value of the first buffer appears on the display.
 - Place the electrodes in the first buffer.
 - Press short 'MODE/RUN'.
 - On the display the pH value of the second buffer appear.
 - Place the electrodes in the second buffer.
 - Press short 'MODE/RUN'.
 - On the display appears CAL^{pH}----mV.
 - Wait till on the display the slope of the calibration curve appears and write it down in the appropriate list.
 - Press short 'MODE/RUN'.
 - The asymmetric potential of the calibration curve appears on the display.
 - Write it down in the appropriate graph.
 - Press long on 'MODE/RUN'.
- On the display 'pH' appears.
- The pH measurements can start now.
- Measure a pH 2 standard and check the value.

3.2.3. pH measurement

- Place the electrodes in the sample.

- Press 'AUTOREAD'.

On the display 'AUTO.R' appears.

- Press shortly 'MODE/RUN'.

On the display ----- appears. and when the value is stable. the pH value of the sample will appear.

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1281 x

1863 apr. '88 A. v.d. Toorn/
ir. K. Nugroho/
drs. J. Harmsen

Methods for the determination of chemical
composition of water samples from column
experiments
(English and Indonesian version)

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4.1. WTW laboratory meter

4.1.1. Calibration

- Take the 0.05 m KCl solution.
- Measure the temperature of the solution.
Adjust the temperature on the meter.
- Rinse the cell with some KCl solution.
- Measure the conductivity. Adjust with CONST ADJ to read 6.66.
- Note the value of CONST ADJ in the control table.
- Measure the conductivity of the 0.01 m KCl solution.
The reading must be 1.41 ± 0.05 .
If not, make new solutions.
Note the value in the control table.

4.1.2. Measurement of samples

- Rinse the electrode.
- Measure the temperature of the samples.
- Adjust this value on the instrument with TEMP COMP.
- Measure the conductivity of the samples. Note the values and also the name of the instrument.
There is no need to correct the values after temperature adjustment.
This is done automatically by the instrument.

4.2. Field equipment

4.2.1. Calibration

- Take the 0.05 m KCl solution.
- Measure the temperature of the solution.
- Rinse the cell with some KCl solution.
- Measure the conductivity. Adjust with CELL to read the value presented in Table 2.
- Measure the conductivity of the 0.01 m KCl solution.

The reading must be the value given in Table 2 with deviation of ± 0.05 .

If not make new solutions.

Note the value in the control table.

Example of Control table for Analytical Control instrument

name	date	const adj	0.01 m KCl	remarks
Antoni	6/4/88	1.10	1.39	new inst new sol

Table 2. Specific conductance of KCl solutions (mS/cm=mho/cm) as a function of temperature (°C)

temperature (°C)	0.01 m KCl	0.05 m KCl
15	1.147	5.404
16	1.173	5.527
17	1.199	5.651
18	1.225	5.775
19	1.251	5.889
20	1.278	6.024
21	1.305	6.149
22	1.332	6.275
23	1.359	6.402
24	1.386	6.529
25	1.413	6.656
26	1.440	6.784
27	1.468	6.912
28	1.496	7.041
29	1.524	7.170
30	1.552	7.300
31	1.580	7.430
32	1.609	7.561
33	1.637	7.692
34	1.666	7.824
35	1.695	7.956

4.2.2. Field equipment measurement

- Rinse the electrode.
- Measure the temperature of the samples and note this on the report.
- Measure the conductivity of the samples. Note the values and also the name of the instrument.

- Correct the values for the temperature as follows:

Conductivity (mS/cm 25°C) = factor x reading.

The correction factors for the temperature is given in Table 3.

Table 3. Temperature correction factors for conductivity

Temperature (°C)	Factor	Temperature (°C)	Factor
15	1.246	26	0.979
16	1.211	27	0.958
17	1.189	28	0.938
18	1.161	29	0.919
19	1.135	30	0.899
20	1.110	31	0.88
21	1.084	32	0.86
22	1.063	33	0.84
23	1.041	34	0.81
24	1.020		
25	1.000		

Use the Analytical Control instrument results as a reference. Every month measure on this meter as well as on all other meters 5 samples. Calculate the relative difference (%) as (reading meter/reading reference -1)x100%. Take the average of the 5 samples. Note this result in a control table. The absolute difference must be smaller then 3%. Unless cleaning of the electrode with HCl can improve the results.

4.3. Chemicals

KCl solution 0.05 m:

KCl solution 0.01 m.

5. DETERMINATION OF CLORIDE CONTENT USING ELECTRODE

5.1. Water samples

- Pipette 20 ml sample in a beaker of 100 ml.
- Add 1.25 ml 1 n HNO_3 .
- Add a magnetic stirrer.
- Dilute to ca 60 ml with distilled water.
- Titrate with 0.05 n AgNO_3 till a potential of -120 mV is read on the mV-meter.
- Note the amount of AgNO_3 used (a ml).
- Note the normality of the solution.

5.2. Calculation

$$[\text{Cl}^-] = 50 \times n \times a \text{ (meq/l)} = [\text{Cl}^-] = a \times n \times 35.5 \times 50 \text{ (mg/l)}$$

n = normality of AgNO_3 solution

a = amount of AgNO_3 used for sample

If another sample volume then 20 ml is used. multiply with; 20/sample volume (ml).

5.3. Chemicals

- 1.0 n HNO_3 : Dilute 7 ml HNO_3 (65%) to 1 liter
- 0.05 n AgNO_3 : Dissolve 8.4935 g AgNO_3 in 1 liter of distilled water (measuring flask).

6. DETERMINATION OF BICARBONATE AND CARBONATE CONTENT

6.1. Water samples

- Pipette 50 ml sample in a beaker of 100 ml.
- Add a magnetic stirrer.
- Measure the pH with a calibrated pH-meter and note the value on the report.
- If the pH is more the 8.3, titrate with 0.1 HCL till this pH is reached.
- Note the amount of HCl used (a ml).
- Titrate with 0.1 n HCl till a pH of 4.1 is reached.
- Note also this amount used (b ml).
- Note the normality of the solution.

6.2. Calculation

$$\begin{aligned} -[\text{CO}_3^-] &= 40 \times n \times a \quad (\text{meq/l}) \quad \text{or} \quad 40 \times n \times a \times 30 \quad (\text{mg/l}) \\ -[\text{HCO}_3^-] &= 20 \times n \times (b-2a) \quad (\text{meq/l}) \quad \text{or} \quad 20 \times n \times (b-2a) \times 61 \quad (\text{mg/l}) \end{aligned}$$

n = normality of HCl solution

a = amount of HCl used to reach pH 8.3

b = amount of HCL used to reach pH 4.1

If another sample volume then 50 ml is used, multiply with; 50/sample volume (ml).

6.3. Chemicals

0,1 n HCl: dilute 8.3 ml of concentrated HCl till 1 liter with distilled water (measuring flask).

7. PHOTOMETRIC DETERMINATION ON NITRATE CONTENT

Salicylic acid is nitrated in a sulphuric acid medium by nitrate. The produced nitro-compound has a intensive yellow colour in a alkalic medium. The absorbance at 415 nm is a measure for the amount of nitrate.

7.1. Apparatus

A spectrophotometer suiteable for measurements at 415 nm. Measuring cells with a pathway of 10 mm.

7.2. Sample

Samples for nitrate can be stored on a cool place (4°C) for one week. If biological activity will occur, the sample must be acidified with H_2SO_4 till pH 2.

7.3. Procedure

- Pipet v ml of sample containing not more than 0.5 mg NO_3 in a beaker of 50 ml.
- Add 1 drop of fenolftaleine and add as much as 0.1 M NaOH till the red colour appears.
- Add 2 ml of sodium-salicylaat reagents.
- Dry on a waterbath (100°C) till the sample is completely dry.
- Add 2 ml of sulphuric acid (18 mol/l) and be sure that all the salt crystals are wetted by the H_2SO_4 .
- Leave it for 10 minutes.
- Add some water and transfer it to a volumetric flask of 100 ml.
(on this moment the solution must be clear. If there is a yellow colour it is necessary to do a blanc of the sample to correct for this colour. This sample blanc has to contain the same amount of sample and all the reagents except sodium-salicylate. The absorbance must be substrated from Ax.).

- Dilute till 70 ml and add carefully 15 ml 10 M NaOH and mix.
- Leave it to cool and dilute till 100 ml.
- After 30 minutes measure the absorbance (Ax) at a wavelenght of 415 μ m against distilled water.
- Measure also the absorbance of the blanc (A bl).
- Calculate the corrected absorbance $A_x = A_x - A_{bl}$.
- Pipette from the stock nitrate solution 0, 1, 2, 3, 4 and 5 ml in a beaker of 50 ml. Treat this as a sample and measure the absorbance against distilled water. Make a calibration curve with the amount of ml's on the horizontal and the corrected absorbance on the vertical axes. The curve has to be lineair.

7.4. Calculations

- Read from the calibration curve the slope (absorbance/ml).
- Calculate the calibration factor $f = \frac{20}{\text{slope}}$
- Calculate the concentration according:

$$[\text{NO}_3\text{-N}] = \frac{f \cdot A_x}{V_o}$$

$[\text{NO}_3\text{-N}]$ = nitrate concentration (mg/L $\text{NO}_3\text{-N}$)

V_o = volume sample in ml's

f = calibration factor

A_x = corrected absorbance

$[\text{NO}_3]$ = 4.429 p $\text{NO}_3\text{-N}$

7.5. Chemicals

- Sulphuric acid 18 mol/L (p=1.84 g/L)
- Sodium-hydroxide 10 mol/L. Dissolve carefully 400 g NaOH in 600 ml water. Dilute till 1000 ml.
- Sodium-salicylatesolution. Dissolve 0.5 g sodium-salicylate in 100 ml water. This solution must be prepared dayly.

- Stock nitrate solution 0.1 mg/ml $\text{NO}_3\text{-N}$. Dry potassium-nitrate KNO_3 for 1 hour at 120°C . Weight 722 mg and dissolve in water in a measuring flask (1000 ml).
- Standard nitrate solution. Dilute 20 ml of the stock-solution till 100 ml in a measuring flask. This solution contains 0.02 mg/ml $\text{NO}_3\text{-N}$.

7.6. References

NEN 6440 (1981). Water photometric determination of nitrate content.
Nederlands Normalisatie Instituut Delft.

8. PHOTOMETRIC DETERMINATION OF SULPHATE CONTENT

Sulphate is precipitated in hydrochloric acid medium with bariumchloride to form bariumsulphate. This is kept in suspension with a conditioning reagent. The absorbance of light by the suspension is measured with a spectrophotometer. The sulphate concentration is established by comparison of the reading with those of standard solutions (standard curve).

8.1. Instrumentation

- Spectrophotometer for use at a wavelength of 420 nm. Cell length 1 cm.
- Magnetic stirrer.
- Stopwatch.
- Erlenmeyer flasks of 100 ml.

8.2. Reagents

- Conditioning reagent.
- Barium-chloride crystals.
- Standard sulphate solution 1000 mg/l.

8.3. Procedure

This procedure is valid for the standard solutions as well as for the samples.

- Pipette 20 ml of the sample in a erlenmeyer flask of 100 ml.
- Add 2 ml of conditioning reagent.
- Add a magnetic rod and place the flask on a stirrer. Place a funnel on the flask.
- Add 0.10 g barium-chloride.
- Stir exactly during 1 minute. Adjust during this period the spectrophotometer on 0.000 using distilled water.

- Measure the absorbance on the spectrophotometer at 420 nm. Read the during 4 minutes. Take the highest reading obtained during this period. If the absorbance is still getting higher, take the absorbance after 4 minutes (this is 5 minutes after adding the barium-chloride). If the reading is more then the reading of the highest standard, dilute the sample.
- Measure a blanc, using 20 ml of the sample, add all the reagents except barium-chloride. Read the absorbance at 420 nm. Depending on the color of the sample the reading will be 0.000 or more.

8.4. Calculation

Calculate the absorbance caused by sulphate;

$$A = A_s - A_b$$

A = Absorption caused by sulphate

A_s = Absorption of sample with BaSO_4

A_b = absorption without BaSO_4 (blanc)

Make a calibration curve as shown (figure 4).

Absorbance


Fig. 4. Calibration curve of sulphate

Read the amount of sulphate from the curve. Multiply it with the dilution used.

8.5. Chemicals/reagents

- Conditioning reagent.

Dissolve in 800 ml distilled water: 75 g NaCl. Add 30 ml concentrated hydrochloric acid. 100 ml 95% ethanol and 50 ml glycerine (glycerol) mix well.

- Barium-chloride crystals.

- Standard sulphate solution 1000 mg/l. Weight 1.479 Na_2SO_4 and dissolve in 1 liter water (measuring flask). The Na_2SO_4 must be dried at 300°C for 4 hours. Pipette from this solution 1, 2, 3, 4, 5, 6, 7, 8, 9 and 10 ml in measuring flasks of 100 ml. This gives solutions of 10, 20, 30, 40, 50, 60, 70, 80, 90 and 100 mg/l.

9. DETERMINATION OF SODIUM AND POTASSIUM WITH FLAMEPHOTOMETER

9.1. Calibration of the meter (monthly)

- Prepare standard solutions with 10, 20, 30, 40 and 50 mg Na^+ and K^+ per litre. For this purpose, dilute 5, 10, 20, 30, 40 and 50 ml of the 1000 mg $\text{NaCl} + \text{KCl}$ solution to 1 liter with distilled water.
- Turn on the meter and light the flame. The flame has to be blue.
- Install the sodium or potassium filter.
- Allow a warm-up period.
- Aspirate alternatively distilled water and the 50 mg/l solution. Operate the sensitivity controls untill the reading of the distilled water is zero and the reading of the standard solution is 100.
- Aspirate all the standard solution and note the readings.
- Make a calibration curve. X-axis = concentration, Y-axis = reading.

9.2. Measurements

- Use the standard solutions with 20 and 50 mg/l.
- Turn on the meter and light the flame. The flame has to be blue.
- Install the sodium or potassium filter.
- Allow a warm-up period of 10 minutes.

9.3. Aspirate

- Alternativily distilled water and the 50 mg/l solution and operate the sensitivity controls untill the reading of the distilled water is zero and the reading of the standard solution is 100.
- Aspirate the 20 mg/l solution. Compare the reading of the standard with the value in the calibration curve. A deviation of 3% is allowed. If the deviation is more then 3%, check the solutions used. If no mistake is found, prepare new solutions. If necessary a new calibration curve.

- Aspirate the samples. Note the reading. If the reading is more than 100, dilute the sample 10 times. Note the dilution which is made.
- Check the performance of the photometer every 5 samples by aspirating the blank and the 50 mg/l solution. Adjust the sensitivity as necessary. If the reading of the highest standard has changed considerably inspect the nebuliser and tubing for clogging and clean it if necessary.

9.4. Calculation

- Read the concentration from the calibration curve. It is also allowed to use a mathematical function describing this curve.
- Multiply the value from this curve with the dilutionfactor.
- If only low concentrations < 10 mg/l are measured, make a calibration curve in this concentration region. An advantage of this low range is, that the calibration curve is almost linear.
- If a mathematical description of the calibration curve is used, be sure that this function is also changed when the calibration curve is changed.

9.5. Chemicals

- Stock solution of 1000 mg Na⁺ and 1000 mg K⁺ per liter.
 - Solute 2106.5 mg NaCl and 1782.6 mg KCl in a measuring flask of 1 liter and dilute to the mark.
- NaCl and KCl have to be dried at 105°C during 3 hours.

10. DETERMINATION OF NA, K, CA, MG, FE AND MN WITH THE ATOMIC ABSORPTION SPECTROPHOTOMETER

1. Empty the bottle on the ground.

2. Set the lamp that you want to use in the right position.

If there have to be another lamp in the lamp-holder, loose the lamp-holder with lamp, by pulling it to the left.


Remove the lamp by pressing the white knob at the backside of the lamp-holder and pull out the lamp.

Place the desirable lamp in the holder by pressing the white knob at the back of the lamp-holder and push the new lamp in the lamp-holder.

Replace the lamp-holder and lamp.

3. Turn on the instrument, 3 knobs power.

4. Press from 'operating lamp' (red knobs) the knob with the number of the lamp that you want to use (see figure 5).

NOTE

1. FIRST POWER SWITCH
2. SECOND POWER SWITCH
3. THIRD POWER SWITCH
4. FINE ADJUSTMENT SWITCH
5. CURRENT ADJUSTMENT DISPLAY
6. CURRENT ADJUSTING SWITCH
7. LAMP STATIC
8. SWITCH TO OPEN THE LAMP CLAMP
9. LAMP NUMBER INDICATOR
10. BURNING ADJUSTMENT
11. ASPIRATOR
12. NEBULIZER
13. BLIT

14. WAVE LENGTH DISPLAY
15. MONOCHROMATOR
16. A-SET. FINE ADJUSTMENT
17. A-SCAN.
18. ACETYLENE AND AIR SUPPLY
19. ACETYLENE AND AIR SYSTEM
20. ACETYLENE AND AIR ADJUSTMENT
21. ACETYLENE AND AIR AMOUNT USED
22. ACETYLENE AND AIR DISPLAY
23. 'READ' AND 'INT' DISPLAY
24. READING DISPLAY
25. READ SELECTION SWITCH
26. EXHAUST SWITCH

Fig. 5. Atomic Absorption Spectrophotometer scheme

5. Set lamp current switch on 'OPR' and adjust with 'FINE OPR' switch the desired current (see table 4).
6. Exhaust on 'MAX'.
7. Adjust the slit (left of the monochromator) (see table 4).
8. Adjust wavelength with 'A scan' (coarse) and 'A set' (fine).
9. Press knobs 'Abs' and 'Damp A'.
10. With 'A set' adjust the wave-length precisely (minimum absorbance).
If necessary adjust the reading to 0.000 absorbance by pressing knob 'zero' till the light goes out.
11. Adjust the lamp to minimum signal with the knobs left on the lamp-holder.
12. If necessary readjust the wave-length with 'A set'.
13. Open the acetylene cylinder (pressure 0.7 atm).
14. Turn on the compressor and adjust the pressure (pressure 3 atm).
15. Support on 'Air' ('Support Press' on 25).
16. Fuel on 'A cet' and ignite the flame.
17. Fuel on desired value (see table 4).
18. Aspirate distilled water.
19. Wait for 15 minutes.
20. Aspirate a blanc and set on zero by pressing 'ZERO' till the light goes out.
21. Aspirate the alignement solution.
Adjust the burner (maximum absorbance).
Right knob, turning the burner.
Left knob, burner forewards or backwards.
Middle knob, height of the burner.
22. Adjust the flame. Fuel lowest knob (not more than 3.3 and not less than 1.5) (maximum absorbance).
23. If necessary readjust the burner.
24. Press knob 'INT' 3. The instrument is intergrating the signal then for 3 seconds.
25. Set blank on zero.
26. Measure the standard series.
For each measurement press knob 'READ'. After the desireable intergration time the redding appears at the display.
27. Measure the samples. After every 4 samples measure the blank and adjust if necessary.

28. Measure after 4 samples also a standard. The reading must be the same as the previous reading (within 3%). If not check on clogging and start again.
29. After the last sample, aspirate distilled water for about 2 minutes.
30. Turn off the instrument.
 1. 'FINE OPR' to the left.
 2. 'POWER' off 3 times.
 3. 'FUEL' off and wait till the flame goes out.
 4. 'SUPPORT' off. Turn of the compressor.
 5. If the instrument is not used anymore that day, close the acetyline cylinder and cover the instrument with a dustcover.

Table 4. Wave-length, slit width, lamp current, acetyline pressure, working area, for the measurement of Na, K, Ca, Mg, Fe and Mn with the atomic absorbance spectrophotometer

element	Na	K	Na	Mg	Fe	Mn
wavelength (n.m.)	589.0	766.5	422.7	285.2	248.3	279.5
slit width	50	300	150	25	50	50
lampcurrent (mA)	10	10	5	3	5	5
acetylene (ml/min)	1.8	2.5	2.2	1.8	2.0	2.3
working area (mg/l)	0-2	0-2	0-10	0-4	0-20	0-5
alignement solution (mg/l)	2	2	10	4	10	5
addition	Ce	Ce	La	La	HNO ₃	HNO ₃ . (65%)

Remarks: Normal working area for Mg is 0-0.4 mg/l. By turning the flame it is possible to work in the range 0-4 p.p.m. Turn the burner head on such a way that the alignment solution gives a reading of about 0.8 absorbance.

Additions:

Ce = 2 ml cesium solution in 100 ml diluting water

cesium solution = 63.4 gr CeCl m 1 liter distilled water

La = 5 ml lanthanum solution in 100 ml diluting water

lanthanum solution = 115 gr LaCl₃ m 1 liter in HCl

SrCl₂ can also be used (180 g instead of 115 g lanthanum chloride).

HNO₃ = 1 ml 65% HNO₃ m 100 ml diluting water.

If the sample has not to be diluted, add one drop 65% HNO₃ in each sample tube from 10 ml.

11. PHOTOMETRIC DETERMINATION OF ALUMINIUM CONTENT

Aluminium produces with pyrocatechol-violet a red/violet coloured solution. The absorbance of the solution measured at 580 nm is a measure for the amount of aluminium.

11.1. Apparatus

A spectrophotometer suitable for measurements at 580 nm. Measuring cells with a pathway of 10 mm.

11.2. Samples

Samples for the determination of aluminium must be clear. If not clear it has to be filtrated. Use in that case a membrane filter 0.45 μ m. Do not use the first part of the filtrate. The noncontent is not allowed to exceed 70 μ g m in the pipetted amount of sample.

11.3. Procedure

- Pipet V_0 ml of sample in a polythene measuring flask of 50 ml.
- Add $0.01 \times V_0$ ml HCl 5 m.
- Add one by one:
 - 1 ml 1.10 fenantroline solution;
 - 2 ml pyrocatecholviolet solution;
 - 10 ml hexamethylenetetramine solution.
- Mix ater every addition. Dilute till 50 ml (the mark) and mix three times.
- Set the wavelength of the spectrophotometer at 580 nm and adjust on zero with distilled water.
- Be sure that the pH is 6.2.
- Measure between 30 and 40 minutes after the addition of the reagents, the absorbance at 580 nm.

- Measure a blank, using V_0 ml distilled water and treat this as a sample.
- Calculate the corrected absorbance using:

$$A_x' = A_x - A_0$$

A_x = absorbance of the sample

A_0 = absorbance measured in distilled water.

Note: If the sample is coloured or it is not quite clear, measure the absorbance A_y , on the same way, as the sample, without using pyrocatecholviolet solution.

The final absorbance is:

$$A_x' = A_x - A_y - A_0$$

- Prepare simultaneously a standard serie of 0, 5, 10, 25 μg aluminium in the analysesample.
- Take for this 0, 5, 10 and 25 ml from the standardsolution II and draw up a calibration curve.

11.4. Calculation

Calculate the aluminium content with:

$$[\text{Al}] = \frac{1000 \cdot b}{V_0}$$

In which: $[\text{Al}]$ = the aluminium content in $\mu\text{g/l}$

b = the amount Al, which is corresponding with the corrected A_x according to the calibration curve

V_0 = the amount of water in ml

11.4. Chemicals

- Hydrochloric acid 5 mol/l.
- Hydrochloric acid 0.1 mol/l.
- 1.10 fenantroline/hydroxylammonium solution. Solute 50 ± 0.5 g hydroxyl-ammonium in about 400 ml water. Add to the solution 0.5 ± 0.005 g 1.10 fenantrolinehydrochlorine. Solute in a measuring flask of 1 liter, fill it to the mark and mix. The solution is allowed to be light/pink. The solution is fenable for two months in a polythene fask.
- Chemicals.
Haxamethylene buffer solution.
Solute 150 ± 0.5 g hexamethiline-tetramine in about 350 ml water. If there are crystals left behind, filter it through fibreglass filtering paper (Whatman GF/C or comparable quality). Bring the solution in a measuring flask of 500 ml and add 6 ml 25% ammonia. The solution is fenable for two months, provided that it is kept in a polythene flask. After use close it immediately.
- Aluminium stock solution 0.5 g/l.
Solute 8.79 g potassium aluminium sulphate ($\text{KAl}(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$) in distilled water. Add 20 ml hydrochloric acid (12 mol/l). Put this in a measuring flask of 1 liter. Dilute with distilled water to the mark and mix.
- Aluminium standard solution I (5 mg/l).
Pipette 10 ml from the aluminium stock solution in a measuring flask of 1 liter and add 20 ml hydrochloric acid (5 mol/l), dilute with distilled water to the mark and mix.
- Aluminium standard solution II (0.1 mg/l).
Pipette 10 ml of the aluminium standard solution I (5 mg/l) in a measuring flask of 500 ml. Add 10 ml hydrochloric acid (5 mol/l) and dilute to the mark with distilled water and mix.
- Both aluminium standard solutions have to be made fresh every day.

12. ADJUSTING NORMALITIES

12.1. Normality of AgNO_3 12.1.1. Normality of 0.005 n AgNO_3

- Pipette 1 ml of 0.05 n NaCl standard solution.

Treat this as a normal sample. Note the amount necessary to reach the end point.

The normality is:
$$\frac{0.05}{\text{Amount titrated (ml)}}$$

12.1.2. Normality of 0.05 n AgNO_3

- Pipette 10 ml of 0.05 n NaCl standard solution.

Treat this as a normal sample. Note the amount necessary to reach the end point.

The normality is:
$$\frac{0.05}{\text{Amount titrated (ml)}}$$

12.2. Normality of HCl

12.2.1. Normality of 0.01 n HCl

- Pipette 2 ml of 0.05 n Na_2CO_3 standard solution.

Treat this as a normal sample. Note the amount necessary to reach the endpoint.

The normality is:
$$\frac{0.1}{\text{Amount titrated (ml)}}$$

or

- Titrate 20.00 ml HCl in a beaker of 100 ml.

Titrate this solution with 0.05 n AgNO_3 to measure the chloride.

Note the volume necessary to reach the end point.

The normality is:

$$\frac{\text{Volume AgNO}_3 \times \text{normality AgNO}_3}{20.00}$$

12.2.2. Normality of 0.1 n HCl

- Pipette 20 ml of 0.05 n Na_2CO_3 standard solution.

Treat this as a normal sample. Note the amount necessary to reach the endpoint.

The normality is: $\frac{1.0}{\text{Amount titrated (ml)}}$

or

- Titrate 5.00 ml HCl in a beaker of 100 ml.

Titrate this solution with 0.05 n AgNO_3 to measure the chloride.

Note the volume necessary to reach the endpoint.

The normality is:

$$\frac{\text{Volume AgNO}_3 \times \text{normality AgNO}_3}{5.00}$$

13. CHLORIDE DETERMINATION USING INDICATOR

13.1. Watersamples

- Pipette 20 ml sample in an erlenmeyer flask of 100 ml.
- Add 10 drops (about 0.2 ml) indicator solution.
- Titrate with 0.05 n AgNO_3 to a pinkish/yellow endpoint.
- Note the amount of AgNO_3 used (a ml).
- Do a blank titration using 20 ml distilled water.
- Note also this amount used (b ml).
- Note the normality of the solution.

13.2. Calculation

- $[\text{Cl}] = 50 \times n \times 35.4 \times (a-b) \text{ (mg/l)}$

n = normality of AgNO_3 solution

a = amount of AgNO_3 used for sample

b = amount of AgNO_3 used for blank

- If another sample volume than 20 ml is used, multiply with;
20/sample volume (ml).

13.3. Chemicals

- Indicator potassium-chromate: dissolve 5 g potassium-chromate in 50 ml distilled water.

14. BICARBONATE AND CARBONATE DETERMINATION USING TITRATION METHODS WITH INDICATOR

14.1. Watersamples

- Pipet 50 ml sample in an erlenmeyer flask of 100 ml.
- Add 1 drop phenolphateline indicator solution.
- If the solution has a red color titrate with 0.1 n HCl till the color dissapears.
- Note the amount of HCl used (a ml).
- Add 2 drops methyl-orange indicator solution.
- Titrate with 0.1 n HCl till the color changes from yellow to orange.
- Note also this amount used (b ml).
- Nota the normality of the solution.

14.2. Calculation

- $[\text{CO}_3^{2-}] = 40 \times n \times a \text{ (meq/l)}$
 $= 40 \times n \times a \times 30 \text{ (mg/l)}$
- $[\text{HCO}_3^{-}] = 20 \times n \text{ (b-2a) (meq/l)}$
 $= 20 \times n \text{ (b-2a) } \times 61 \text{ (mg/l)}$

n = normality of HCl solution

a = amount of HCl used till changing red color

b = amount of HCl used till changing yellow color

14.3. Chemicals

- Phenolphateline
- Methyl-red

15. CALCIUM AND MAGNESIUM DETERMINATION USING TITRATION

This method is suitable to determine calcium and magnesium with a minimum concentration of 1 mg/l.

15.1. Calcium and magnesium (total hardness)

15.1.1. Procedure

Pipet W_1 ml in a conical flask and dilute with distilled water to ± 100 ml. Add one by one : 5 ml buffer solution; 5 ml triethanolamine solution; 5 ml magnesium EDTA solution; 10 drops (± 0.5 ml) potassium-cyanide and as much as indicator that the solution is coloured. Mix after every addition. Titrate with the Na_2 EDTA solution, till the green colour just does not change anymore (use V_1 ml).

15.1.2. Calculation

$$H_t = \frac{1000 V_1 c}{W_1}$$

H_t = total hardness in mmol/liter

V_1 = used amount of Na_2 EDTA (ml)

c = molarity Na_2 EDTA solution (mol/liter)

W_1 = pipetted amount of water used with calcium determination

15.2. Calcium determination

15.2.1. Procedure

Pipette ml (W_1) in a conical flask and dilute with distilled water to ± 100 ml. Add one by one : 5 ml triethanolamine solution; 4 ml potassiumhydroxide solution, ten drops (about 0.5 ml) potassium-cyanidesolution. and as much as indicator that the solution is coloured red. Mix after every addition. Titrate with the Na_2 EDTA solution, till the blue colour just does not change anymore (use V_2 ml).

15.2.2. Calculation

$$C \text{ Ca}^{2+} = \frac{1000 V_2 c \times 40.08}{W_2}$$

In which: $C \text{ Ca}^{2+}$ = content of calcium in mg/l

V_2 = amount of Na_2 EDTA used with calcium determination (ml)

c = molarity Na_2 EDTA solution (mol/l)

W_2 = pipetted amount of water used with calcium determination

15.3. Magnesium determination

Magnesium content is determined after total hardness and calcium content determination.

$$C \text{ Mg}^{2+} = (H_t - \frac{C \text{ Ca}^{2+}}{40.08}) \times 24.32$$

$C \text{ Mg}^{2+}$ = magnesium content in mg/liter

H_t = total hardness

$C \text{ Ca}^{2+}$ = calcium content in mg/liter

Remarks

If the water contents more than 10 meq carbonate and/or hydrogencarbonate a liter, add some nitric acid (till pH 2), soil short, and neutralize on methylred.

If the pipetted sample contents more than 1.5 mg iron, the added triethanolamine solution is not sufficient to neutralize the disturbance. In that case add 5 ml triethanolamine solution for every 2 mg iron.

15.4. Chemicals

- Na_2EDTA solution 0.025 mol/l (molarity c)

Solute 10 g $\text{Na}_2\text{EDTA} \cdot 2\text{H}_2\text{O}$ in a measuring flask and dilute with distilled water till 1 liter. Keep this solution in a polythene bottle.

- Standard calcium solution

Solute 1000 mg (= a) calcium-carbonate in 50 ml hydrochloric acid 5 m, boil the solution for some minutes and dilute after cooling in a measuring flask of 500 ml. Before diluting dry the calcium-carbonate for 2 hours at 400°C.

- Indicators

For hardness determination.

Mix 0.1 g eriochrom-black -T and 0.04 g methylred with 30 g potassium-chloride in a mortar.

For calcium determination.

[2 hydroxy-1-(2 hydroxy-4-sulfo-1-naftylazo)-3 naftalene carbonacid with 30 g potassiumchloride in a mortar.

- Potassium-hydroxide solution

Solute 50 g potassium-hydroxide in 100 ml distilled water.

- Magnesium EDTA solution (3 g/l)

Solute 3 g $\text{MgEDTA} \cdot 6\text{H}_2\text{O}$ in 1 liter distilled water.

- Buffersolution pH 10.3-10.5

Solute 75 g glycine ($\text{H}_2\text{NCH}_2\text{COOH}$) in 500 ml distilled water. Add as much as potassium-hydroxide solution till the pH is 10.3-10.5. Keep this solution in a polythene bottle.

- Potassium-cyanide solution (toxic)

Solute 10 g potassium-cyanide in 100 ml distilled water.

- Titration Na_2 EDTA solution 0.025 mol/l

1 Pipet 50 ml of the calcium solution in a conical flask of 250 ml. dilute with distilled water to 100 ml, add 5 ml buffer solution, 5 ml triethanol solution, 5 ml Mg.EDTA solution, 10 drops (about 0.5 ml) potassiumcyanide solution, and as much as indicator that the solution is coloured. Titrate slowly with the Na_2 EDTA solution till the green colour just is not changing anymore.

Molarity $c =$

$$c = \frac{a}{10 M_{\text{CaCO}_2} \cdot V}$$

in which: M_{CaCO_2} = molar mass of calciumcarbonate (100.1 g/mol)

c = molarity of Na_2 EDTA solution

a = weight amount of calciumcarbonate (mg)

V = volume used of Na_2 EDTA solution (ml)

- Pipet 50 ml of the calcium solution in a conical flask of 250 ml. dilute till 100 ml with distilled water, add 5 ml triethanolamine-solution, 4 ml potassium-hydroxide solution, 10 drops (about 0.5 ml) potassium-cyanide solution and as much as indicator that the solution is coloured. Titrate with Na_2 EDTA solution, till the green colour just is not changing anymore.

Calculation

$$c_1 = \frac{a}{10 M_{\text{CaCO}_2} \cdot V}$$

in which: M_{CaCO_2} = molar mass of calciumcarbonate (100.1 g/mol)

c_1 = molarity of Na_2 EDTA solution (mol/l)

a = weight amount of calciumcarbonate (mg)

V = volume used of Na_2 EDTA solution (ml)

c and c_1 must give results within the limits of duplicate measurements.
If not check the chemicals used.

Chemicals:

Phenolphthaline indicator: dissolve 1 g phenolphthaline in 100 ml ethanol 96%.

Methylorange indicator: dissolve 100 mg methylorange in 1 liter distilled water.

0.1 n HCl Dilute 8.3 ml of concentrated HCl till 1 liter with distilled water (measuring flask).

LITERATURE

AMERICAN PUBLISHING. Health association standard methods for the examination of water and waste water 14th edition. American Publ. Health Assn. Washington D.C.

NEN, 1979. Determination of total hardness and calcium and magnesium content NEN 6441.

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