Afdeling Algemene Chemie 1986-01-20

Rapport 86.12 Pr.nr. 505.2090

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1981.05.10

AN 30/81

DETERMINATION OF KJELDAHL NITROGEN CONTENT WITH KJELTEC AUTO SYSTEMS I, II, III, and IV

I. INTRODUCTION

The Kjeldahl Nitrogen method, published nearly a century ago, is today the most frequently used procedure when measuring protein content in biological materials. It has also found its application as a method for determination of Nitrogen in many inorganic materials. The reasons for the wide acceptance and its official status are the high reliability in combination with the vast knowledge that has been developed during the long period of use. The Kjeldahl method is also highly appreciated for its general applicability on all types of materials. Originally the method was very tedious and time consuming to perform.

The KJELTEC AUTO SYSTEMS have been developed to meet today's requirements for efficiency, speed and convenience and to reduce the need for chemicals, energy and time, while still allow-

ing the performance of most of the modifications of the original method which are described in literature. This Application Note will in a broad sense describe how determinations can be performed using the different Tecator Digestion Systems in combination with the Kjeltec

Auto 1030 Analyzer. Methods for different products are presented in specific Application Sub Notes.

II. MATERIAL

Dry and well-mixed samples such as seeds, chopped straw, etc. should be ground to a particle size that allows a representative sample of about 0.2 - 1 gram to be taken. Semisolid samples such as meat and meat products, canned or frozen foods, should be homogenized and mixed carefully to allow a representative sample of typically 0.5 - 2 gram to be taken. Solutions, such as fresh milk, may be analysed without pretreatments. Homogenization may, however, be needed if e.g. the cream has separated from milk.

Note: The lower weights are recommended for semimicro Kjeldahl analyses using the Digestion Systems 12 and 40.

In addition the following equipment and accessories are needed for the different Kjeltec Auto Systems:

K J E L T E C A U T O S Y S T E M S				
	I (for macro Kjeldahl)	II	III (for semi-micro Kjeldahl)	IV
1) Digestor with built-in thermostat	Dig.Syst.6 1007-007 (220V/50Hz)	Dig.Syst.20 1015-001 (220V/50Hz)	Dig.Syst.40 1016-001 (220V/50Hz)	Dig.Syst.12 1009-001 (220V/50Hz)
		1007-008 (110V/60Hz)		1009-002 (110V/60Hz)
2) Stand for tubes	1000-0831 (6 tubes)	1000-0082 (20 tubes)	1000-0072 (40 tubes)	1000-0537 (12 tubes)
3) Set of digestion tubes	1000-0155 (6 tubes)	1000-0158 (20 tubes)	1000-0735 (40 tubes)	1000-0825 (12 tubes)
4) Heat shield	1000-0536	1000-0188	1000-0188	1000-0536
5) Exhaust Manifold with water aspirator (incl. stand and drip pan)	1000-0835 (incl. in 1007-007)	1005-034	1006-025	1000-0836
6) Stand for rapid cooling of tubes		1000-0424	1000-0424	
7) Retainer plate		1000-0073	1000-0095	1000-0538
8) Scrubber Unit (optional)		1005-068	1005-068	

- 9) Kjeltec Auto 1030 Analyzer including storage tanks for alkali, receiver solution, separate feed water and titrant (1030-001, 220V/50Hz), (1030-002, 110V/60Hz).
- 10) Special checking thermometer (1000-0103).

III. EQUIPMENT

Suitable equipment for sample preparation is needed. For grains etc. a Cyclotec Sample Mill (1900-0300) or a Cemotec Mill (1090-001) may be used whereas for semisolid samples a Moulinex or Ultra Turrax homogenizer may be used. For weighing of sample the 6002 Analytical Balance is recommended.

6. Digest for 3 - 5 minutes with maximum air flow through the exhaust manifold. Then adjust the flow until fumes are just contained. (If gases tend to escape out into the air even at full air flow, the digestor can be tilted slightly to help avoiding this. Just place a piece of wood about 12 mm underneath each front leg of the digestor. This will prevent condensate from dropping from the collector straight down into the boiling sample mixture, which then may result in sudden and irregular changes in pressure.)
7. Continue digestion until ready. Remove the stand with tubes and exhaust manifold and place the entire assembly in the cooling stand beside the digestor. Increase the air flow through the exhaust manifold if necessary.
8. As soon as the sample solutions have cooled sufficiently they should be diluted with water and mixed (25 ml water for semi-micro and 75 ml for macro tubes). Use heat resistant gloves for protection. (Tap water may be used for this purpose if proven N-free). Cooling can be speeded up by blowing air between the tubes with a small fan. If the digest is too hot when water is added, the reaction will be too violent and sample material might be lost. On the other hand if digest is too cold when water is added, salts may precipitate, which are difficult to redissolve. (Solidification should be avoided but in case precipitation does occur, dissolve it by placing the tubes in the digestor for a short while.)
9. Start up the Kjeltec Auto 1030 Analyzer as described in the manual. Make sure there is no air in the burette or the tubings from the burette.
10. Select Kjeldahl programme and run one or two blanks with the constants set as follows:

$$A = 00.00, \quad B = 1.000, \quad \text{Blank} = .00$$

Note the display value and use it as blank value during the following distillations.

11. Put the prepared* digestion tube in position and close the safety door. Set the B constant to correct value (see table in the manual). B is set to 1.000 if the result should be in ml.
 12. When the "cycle over" lamp is on, note the result and open the safety door (the result is printed if a printer is connected).
- Result (in display) = A + B × (ml titrant - ml blank).
13. Continue from 11 with the next sample.

* If mercury is used as a catalyst and no thiosulphate has been added to the alkali, then 5 ml of the thiosulphate solution (as per 'Kjeltabs M/3.5' used in the digestion) must be added to the sample prior to the distillation.

VI. RESULTS AND COMMENTS

The KJELTEC AUTO SYSTEMS are specially designed to be flexible and thus to cope with most of the Kjeldahl methods used today. The methods presented in this Application Note as well as in the Application Sub Notes for different types of materials are only guidelines showing one or more ways to achieve reliable results. However, the recommendations are based on years of practical experience and collaborative studies together with many of our customers.

In the following we will discuss some of the factors that should be considered when determining nitrogen in organic materials by the Kjeldahl method.

1. Type of sample material

It is well-known that some types of material are far more easy to digest than others. Cereal seeds for example, are known to be simple to digest whereas fish meal and particularly tobacco leaves are far more difficult to digest. Nicotinamide is, because of its resistance to acid digestion, frequently used as a test sample to check the recovery of nitrogen.

2. Sample size

The larger the sample weight is that has to be heated and digested, the more chemicals will be needed, the longer time it will take to reach decomposition temperature of the digest and the longer time will be needed for complete digestion.

Thus it is not recommended to use larger samples than necessary to obtain a representative result.

3. Amount of acid used for digestion

The controlled conditions during digestion eliminates the potentially large loss of acid that may otherwise cause loss of nitrogen. The acid volume used in the KJELTEC AUTO SYSTEMS is therefore generally less than recommended in conventional procedures.

Thus when conventional procedures prescribe the use of 25 to 30 ml of acid it is normally sufficient to use 10 to 12 ml for macro analyses (250 ml tubes) and only about 5 ml for semi-micro analyses (100 ml tubes) when the Kjeltec Auto Systems are used.

4. Choice of catalyst

Extensive research published in literature as well as results from our own labs have shown that mercury is the number one catalyst as long as speed and "complete" recovery are given the highest priority. Its toxicity has stimulated many workers to use other catalysts instead. However, none of them is completely harmless.

In our earlier work we have been using Selenium as the normal catalyst. Kjeltabs with other catalysts are now available. In the Application Sub Notes for the different

materials we sometimes give alternative procedures if the same, or nearly the same, results can be achieved in reasonable time. (If you use Mercury as the catalyst, it can be recovered by collecting the distillation residue from the digestion tubes in a tank where the precipitate can be collected.)

5. Addition of salt to digestion mixture

The effect of salt as an accelerator of the digestion process by increasing the boiling temperature is well documented. Generally Potassium Sulphate is recommended and used because of the positive and long experience of it in Kjeldahl laboratories throughout the world. The acid: salt ratio should, however, be kept under control as too high salt concentrations may cause loss of nitrogen. As a thumb-rule, the salt concentration should not be so high that the digest solidifies when it is cooling down (even without addition of water!!!).

6. Use of oxidizing agents

Hydrogen Peroxide has been and still is used extensively as an oxidizing agent in many laboratories. From time to time, however, small losses of nitrogen may occur, particularly when it is used in combination with Selenium. As a matter of fact we have even found that it is very seldom necessary to use oxidizing agents as a means of reducing the foaming problems. The speed of digestion also seems to be more or less unaffected by oxidizing agents. Therefore we do not recommend the use of Hydrogen Peroxide as standard. However, if it is used, the losses of nitrogen should be carefully checked, especially if the Kjeltabs contain Selenium catalyst.

Hydrogen Peroxide should always be used with care. The violence of the reaction that

takes place when peroxide is added to the sample/acid mixture, depends on

- (a) the composition of the sample
- (b) the concentration of the peroxide used and
- (c) the care taken during the addition.

Particularly when a new batch of Hydrogen Peroxide is to be used or when a completely new type of sample is to be analyzed, special care is recommended.

7. Temperature at which digestion takes place

The digestion temperature is dependent on the boiling point of the digestion mixture, which again is dependent on the acid: salt ratio of the same mixture.

Another factor of significant importance is at what altitude the digestion takes place.

The higher the altitude is, the longer digestion time will be needed due to reduced boiling temperature of the digest. At e.g. 2,000 metres the boiling temperature is about 11°C lower than at the sea level.

The temperature of the digestor is normally set to 420°C, which is well above the boiling point under the analytical conditions recommended. Still it is low enough to avoid the losses of Nitrogen that otherwise may take place if the digests are locally overheated.

8. Duration of Digestion

Depending on the choice of analytical conditions (1) to (7) above, the speed of digestion may vary from approximately 10 minutes up to several hours. If the methods given in the Application Notes for specified type of materials are followed, the digestion will normally be completed within 20 - 45 minutes.