Project 101.6040

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
Section 3

3.2.2.3 DATA Display and Mode Keys

The functions of the Mode keys are explained in Table 3-4. When a mode key is pressed its associated lamp lights, the key function can be cancelled only by pressing another mode key. The first four digits of the display, labelled data, are a numerical value, whose power to the base 10 is indicated by the last two digits, labelled EXP. (Exponent).

DATA	EXP.	VALUE
6.160	BLANK	6.160
6.160	03	$6.16 \times 10^3 = 6160$
6.160	06	$6.16 \times 10^6 = 6160000$

If the DATA display gives the indication U1, U2, U5 or U6 this indicates a defect on the Microcomputer board, which cannot be rectified by the operator. Contact the Service Dept. of Coulter Electronics Ltd., or its appropriate authorized Distributor.

TABLE 3-4 DATA DISPLAY AND MODE KEYS

Key	Data Display on ZM	Printer or VDU Format (80 Characters)
DIA (T _U)	The diameter of the sphere representing the volume to which the upper threshold control is set.	The numerical value to which the upper and lower thresholds are set as well as the diameter of the sphere representing the volumes to which the thresholds are set. If the upper threshold is set to 99.9, no upper threshold data is displayed.
VOL (T _U)	The volume to which the upper threshold is set.	The numerical values to which the upper and lower thresholds are set, as well as the volume of the spheres represented by the threshold settings. If the upper threshold is set to 99.9, no upper threshold data is displayed.
DIA (T _L)	The diameter of the sphere representing the volume to which the lower threshold control is set.	The numerical value of the lower threshold and the diameter of the sphere representing the volume to which the lower threshold is set.
VOL (T _L)	The volume to which the lower threshold is set	The numerical value of the lower threshold and the volumes of the spheres represented by the lower threshold setting.
CAL	Enables the function of the calibration keys and allows new calibration data to be keyed in when the appropriate calibration key is pressed. Overrides Error Code 1, permitting adjustment of aperture current and threshold settings during calibration count, to establish the optimum settings.	When MAN Interface and CAL Modes are selected, all the calibration values are displayed as well as the analysis time (set by the operator, or time taken for a count in manometer mode), the CAL lower threshold numerical setting for the calibration, the aperture current used for calibration and the attenuation switch setting. The sample ID can be entered by using the numerical keys when the CAL mode key is pressed without pressing any of the calibration keys.

3.2.2.4 Data Display and Calibration Keys

When the ZM is not in the CAL mode the operator can display the calibration values that have been keyed in by pressing the appropriate calibration key. The data is then displayed for approximately 1.5 seconds before reverting to its previous mode.

The procedure used for changing the value stored in the ZM for each calibration key is the same. The procedure is as follows:-

- (1) Press the CAL mode key.
- (2) Press the appropriate calibration key, the calibration data stored is then displayed on the data display.
- (3) The new value can then be keyed in using the numerical keys below the data display. If a mistake is made, press the CE key which clears the display only and key in the new value.
- (4) To place the new calibration value in the store, press the appropriate calibration key for the second time.
- (5) When changing one or all of the calibration key values is completed, press any of the Mode keys except CAL to complete the procedure.

The function of the calibration keys is given in Table 3-5.

TABLE 3-5 DATA DISPLAY AND CALIBRATION KEYS

CAUTION: All Data associated with the keys listed in TABLE 3-5, with the exception of TIME(S), are retained in the microprocessor memory when the instrument is switched off. The TIME value is reset to zero to enable manometer control at switch on.

Key	Data Display	Function
DIA (μm)	The assay value of the diameter of the calibration material used during calibration.	To enter the calibration material diameter and to display the value stored in the ZM. The diameter is used to calculate calibration constant.
VOL (fl)	The assay value or calculated value of the calibration material volume.	To enter the calibration material volume (if the calibration material volume is given and not the diameter) and to display the value stored in the ZM. When the key is pressed the diameter is automatically calculated and transferred to the diameter store for calculation of the calibration constant.
MAN VOL (μl)	The manometer volume, selected in the Sampling Stand by the operator, in μl .	To display and enter the volume of electrolyte passing through the aperture during a count when the instrument is operating in the Manometer mode. The manometer volume is required for the ZM calculations of coincidence correction to be applied to the count. If the instrument is operating in the Time Mode, then a calculated value of the volume that passes through the aperture for the selected time must be entered.

Cont'd/...

Section 3

Table 3-5 Continued

Key	Data Display	Function
AP DIA (μ m)	The diameter of the orifice tube aperture used on the Sampling Stand.	To display and enter the diameter of the aperture tube being used in the Sampling Stand. The aperture diameter is required for the coincidence correction calculation.
TIME(S)	If zero is displayed this indicates that the instrument is operating in the manometer mode. If any other figures appear in this display this indicates that the ZM is operating in the Time mode, and the function of the Sampling Stand manometer has been inhibited. Note: In Time mode, maximum time which can be keyed in is 655s.	To display and enter the time period of the count. If zero has been entered, or the timer reset by switching the instrument off, it indicates to the operator and to the ZM that the count is related to a volume of electrolyte as selected by the Manometer Volume Switch on the Sampling Stand. If the ZM is in the Manometer mode and the CAL key is not pressed, the count time is accessible by pressing the TIME key. The count time is then displayed for approximately 1.5 seconds. In the Time mode, the function of this key (when CAL key is not pressed) is to initiate the count when pressed. If the key is pressed during a count, the count stops at the instant it is pressed.
Kd	The diameter calibration factor of the instrument.	When pressed, displays the calibration factor for the previous calibration. The operator can also key in a known factor for a given orifice tube. (once Kd factor is known for an orifice tube it does not have to be recalibrated or Kd re-calculated unless the PRESET GAIN switch setting is changed).

3.2.3 INTERNAL SWITCHES

The internal switches shown in Figure 3-3 are accessible only when the side panels are removed and the front panel door opened. The function of the internal switches is as follows:-


- (1) Security code switch, when set to 'on' inhibits data entry in the CAL mode unless the security code 9876 is first keyed in, then cleared with the CE key.
- (2) The Baud Rate switches are used to make the ZM compatible with a data interface.
- (3) The format switch is used to select either the 80 character format or the 20 character format (See Fig. 2-1 and 2-2).
- (4) The de-limit switches can be used to set any letter in the ASCII code required by the data terminal equipment.

3.2.4 SAMPLING STAND CONTROLS AND INDICATORS (Fig. 3-4)


3.2.4.1 Power Switch

Power OFF, lamp extinguished

Power ON, lamp lit.

THE POSITION OF THE CURSORS
ON THE 'FULL' MENU BECOME
THE LIMITS OF THE 'WINDOW'
MENU


Section 2

2.2 **SAMPLE ANALYSIS**

2.2.1 **FULL RANGE ACCUMULATION**

(1) Press the C256 FULL key the following menu is displayed.


(2) Using the Menu Cursor keys step down the 'Full' menu option and enter the information as follows:

SAMPLE: Enter the Sample Number.

Notes:(1) The previous information is erased directly the first digit is entered.

(2) The flashing cursor remains on each menu in the position it was last placed whilst the instrument is switched on.

= 100%: The number entered here represents the number of particles counted in a channel for 100% display scale. To change the scale press Menu Key $\leftarrow$ to double or $\rightarrow$ to halve. The range available is 125, 250, 500, 1,000, 2,000, 4,000, 8,000, 16,000, 32,000 and 64,000.

ration data
meter
ibration
r or volume
the
a
with the
l settings
r which
displays the
Kd key is
red in the

ZM's memory and is used to calculate and display the size corresponding to any position of the instrument sensitivity controls. Also incorporated is an analysis time control, values for which are entered in the same way as the calibration data, and which can be used to replace the mercury syphon sampling, for very small or very large apertures. When not in the CAL mode these keys are disabled preventing the accidental entry of erroneous data.

ONS

range.
tions to be
eight)

counts per
on at,
et particle
ent on

ing Range
or more.

Timer Range

10ms - 655 seconds

Operating Temp. Range

10°C - 35°C

Operation Humidity Range

10% - 85% Relative Humidity

Data Outputs

Total number of particles above
preset size.
Total number of particles between
preset sizes.
Time of analysis.
Instrument settings.


Interface Output

ASCII, RS232, 20mA current loop.
Variable baud rate 110 - 9600.

npling land	Vac/Waste	C256	MCV/HCT
15.0	11.5	29.0	4.5
17.5	27.0	59.0	34.3
16.0	16.0	39.0	42.0
14.5	2.4	15.0	2.3

.) 220 - 240 V.a.c. $\pm$ 10% 50Hz/60Hz
110 - 120 V.a.c. $\pm$ 10% 50Hz/60Hz

V

The Coulter Channelyzer® 256 is a new, high resolution, size distribution accessory, designed specifically for use with the Coulter Counter® Models ZM and TALL. Size distributions of biological cells and other narrow range material may now be obtained in real time with 64, 128 or 256 channel histograms providing the operator with the highest level of sensitivity in sizing resolution. Size differences as small as 0.05fl may be

detected, making the system an invaluable and necessary tool for everyone involved in areas such as yeast cell growth, red and white blood cell studies for human, veterinary and pharmaceutical applications, cancer research, immunology and some areas of microbiology.

Real Time Analysis

Data in the form of analogue pulses

representing particle number and volumetric size from the Coulter Counter®, are accumulated directly in the memory of the C256 in either 64, 128 or 256 channel resolution as selected by the operator.

COULTER CALIBRATION STANDARD

Assay Sheet

E.S.C.C. CALIBRATOR

BATCH M.10

MATERIAL P.D.V.B. LATEX




DO NOT USE AFTER MARCH, 1989

FOR USE WITH THE COULTER COUNTER®
MODELS MCC, TAIL/PCA, TA, ZM, ZB, ZBI, ZF,
FN, C.1000, P.64/P.128, D. INDUSTRIAL, B & A.

	100 μ m APERTURE	140 μ m APERTURE (MCC)
SINGLET NUMBER MEDIAN DIAMETER (figure A) (For use with MCC, ZM, ZB, ZBI, ZF, FN, D.Ind., B & A)	5.86 μ m	5.85 μ m

Other parameters:

NUMBER MODE (figure B) (For use with C.1000)	5.75 μ m	5.73 μ m
SIZE DISTRIBUTION (ZM/C. 1000) — see overleaf		
NUMBER PEAK SPLIT (figure C) (For use with TA series)	5.80 μ m	5.83 μ m
RECOMMENDED USAGE RATE 100 μ m Aperture, 3–4 drops per 50ml. 140 μ m (MCC) Aperture, 1 drop per 25ml.		

ASSAY VALUES have been determined in SOMATON® diluent prepared by Coulter Electronics Ltd., using COULTER COUNTER® Models TAIL/PCA, ZM and CHANNELYZER® C.1000 which are calibrated and checked with primary reference particles (1).

All particles have a coefficient of variation from their mean size of less than 13% by diameter.

CALIBRATION PARTICLE SUSPENSION: E.S.C.C. (Electronic Somatic Cell Counter) latex calibrator particles are suspended in an aqueous dispersant plus preservative, such that, at the recommended usage rate shown above, the correct concentration (approx. 5% coincidence) for calibration is obtained for direct application to the aperture. Details of the 'cold' (direct) COULTER COUNTER® model M.C.C. calibration technique can be obtained from Coulter Electronics Ltd or their distributors. There is no evidence of E.S.C.C. calibrator latex instability with time, as packed, but for maximum security it is recommended that the product be discarded after the expiry date given above.

SIZE CALIBRATION: The COULTER COUNTER® can be self-calibrated by measuring a suspension of a narrow size range of the material under investigation at a known concentration. From the immersed relative density (g/ml) of the material and from the total volume measured, the calibration factor can be calculated for the COULTER COUNTER® and sensing aperture in use (2). This procedure will remove any doubt concerning the accurate calibration with very irregular or conducting materials.

In practice it is more convenient to calibrate using narrow range particles which have already been sized by a variety of other techniques.

Extensive electron microscopy showed that an original reference latex produced by the Dow Chemical Co. (Midland, Michigan, U.S.A.) with a stated mean size of 3.49 μ m was in fact 3.40 μ m diameter (3,4). Harfield and Wood (1) confirmed this size with a COULTER COUNTER® MODEL B self-calibrated by the integration of data from a sample of alumina, and verified by reference to N.B.S. Standard Reference Material 1033 glass beads, and also by reference to a packed volume of latex particles.

Using aperture tubes and Coulter Counter® models (ZM, C1000, TAIL/PCA) correctly calibrated as above, the E.S.C.C. calibrator is sized by the relevant technique. The single distribution limits for the singlet number median size quoted above (figure A) are defined as one-half and three-halves of the median particle volume. This method of calibration is most convenient for single and dual threshold COULTER COUNTER® models. See instruction manuals for procedure.

A COULTER CHANNELYZER® may be more conveniently calibrated using the modal size of the particles; the 'number mode' (C1000); (figure B).

When calibrating a COULTER COUNTER® model TA, or TAIL without PCA, use the cal. meter method of splitting the number peak; with PCA split the number peak directly (figure C).

10ml

SOMACOUNT®
PRE-FIXED MILK CELL REFERENCE CONTROL FOR THE COULTER COUNTER®

Assay Sheet

BATCH 36 LOW

SEALED VIALS – DO NOT USE AFTER

31st MARCH, 1987

OPEN VIALS – DO NOT USE AFTER

7 DAYS

FOR USE WITH THE COULTER COUNTER®
MODELS M.C.C., ZM, ZF, ZB, ZBI, FN, F, D, IND, B, TAIL, TA & A.

	COULTER® MILK CELL COUNTER (M.C.C.)			OTHER MODELS		
PARAMETER	MEAN VALUE IN SOMATON® (± 1 S.D.)	EXPECTED RANGE A	EXPECTED RANGE B	MEAN VALUE IN SOMATON® (± 1 S.D.)	EXPECTED RANGE A	EXPECTED RANGE B
SOMATIC CELL COUNT* (x10 ⁶ /ml)	0.204 ± 0.011	0.188–0.220	0.176–0.232	0.226 ± 0.011	0.212–0.240	0.203–0.249
MODAL SIZE (µm)	Nominal 5.5µm – size dependent upon analytical method.					

EXPECTED RANGE A: expected range of values of average of 5 samplings from any one vial. (See notes 1 & 2).

EXPECTED RANGE B: expected range of values for any individual measurement. (See note 3).

SIZE VALUE quoted is after 1:100 dilution in SOMATON® and incubation for 10 minutes at 80°C. Expressed as micro-meters (µm) equivalent spherical diameter.

*SOMATIC CELL COUNT ABOVE 4.5µm DIAMETER. (Typical setting when using SOMATON® and SOMAFIX®). To determine assay values for any other diameter (d), as required by the particular fresh milk fixation process in routine use:

TABLE 2

d (µm)		x10 ⁶ /ml	d (µm)		x10 ⁶ /ml
4.0	ADD	0.013	4.6	SUBTRACT	0.002
4.2	ADD	0.007	4.8	SUBTRACT	0.005
4.4	ADD	0.005	5.0	SUBTRACT	0.008

Somacount® can provide a simple visual check of the size calibration. The modal size on the M.C.C. is directly measured in Somaton® at 80°C. Modal size on other models is measured post-incubation in Somaton® at ambient temperature.

Accurate size calibration of the threshold scale is best achieved via the modal or singlet number median sizes of the Coulter E.S.C.C. calibration latex (part number 9900526) using the recommended instrument calibration procedure.

ASSAY VALUES have been determined using COULTER COUNTER® Models M.C.C., ZM and CHANNELYZER® C.1000 which are calibrated and checked with primary reference particles (1). The electrolyte used is SOMATON®, a milk cell diluent prepared by Coulter Electronics Ltd.

FOR IN VITRO DIAGNOSTIC USE ONLY

9942216 3/85

Bijlage 19 b
10ml

SOMACOUNT®

PRE-FIXED MILK CELL REFERENCE CONTROL FOR THE COULTER COUNTER®

Assay Sheet

BATCH 236 HIGH

SEALED VIALS – DO NOT USE AFTER 31st MARCH, 1987

OPEN VIALS – DO NOT USE AFTER 7 DAYS

FOR USE WITH THE COULTER COUNTER®
MODELS M.C.C., ZM, ZF, ZB, ZBI, FN, F, D, IND, B, TAIL, TA & A.

	COULTER® MILK CELL COUNTER (M.C.C.)			OTHER MODELS		
PARAMETER	MEAN VALUE IN SOMATON® (± I.S.D.)	EXPECTED RANGE A	EXPECTED RANGE B	MEAN VALUE IN SOMATON® (± I.S.D.)	EXPECTED RANGE A	EXPECTED RANGE B
SOMATIC CELL COUNT* (x10 ⁶ /ml)	1.375 ± 0.019	1.340-1.410	1.306-1.444	1.499 ± 0.019	1.468-1.530	1.442-1.556
MODAL SIZE (µm)	Nominal 5.5µm - size dependent upon analytical method.					

EXPECTED RANGE A: expected range of values of average of 5 samplings from any one vial. (See notes 1 & 2).

EXPECTED RANGE B: expected range of values for any individual measurement. (See note 3).

SIZE VALUE quoted is after 1:100 dilution in SOMATON® and incubation for 10 minutes at 80°C. Expressed as micro-meters (µm) equivalent spherical diameter.

*SOMATIC CELL COUNT ABOVE 4.5µm DIAMETER. (Typical setting when using SOMATON® and SOMAFIX®). To determine assay values for any other diameter (d), as required by the particular fresh milk fixation process in routine use:

TABLE 2

d (µm)		x10 ⁶ /mL	d (µm)		x10 ⁶ /mL
4.0	ADD	0.008	4.6	SUBTRACT	0.002
4.2	ADD	0.006	4.8	SUBTRACT	0.011
4.4	ADD	0.003	5.0	SUBTRACT	0.021

Somacount® can provide a simple visual check of the size calibration. The modal size on the M.C.C. is directly measured in Somaton® at 80°C. Modal size on other models is measured post-incubation in Somaton® at ambient temperature.

Accurate size calibration of the threshold scale is best achieved via the modal or singlet number median sizes of the Coulter E.S.C.C. calibration latex (part number 9900526) using the recommended instrument calibration procedure.

ASSAY VALUES have been determined using COULTER COUNTER® Models M.C.C., ZM and CHANNELYZER® C.1000 which are calibrated and checked with primary reference particles (1). The electrolyte used is SOMATON®, a milk cell diluent prepared by Coulter Electronics Ltd.

FOR IN VITRO DIAGNOSTIC USE ONLY

9942216 3/85

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Copyright 1990, Rijks-Kwaliteitsinstituut voor land- en tuinbouwprodukten.