



**Project SOP D-JRP17-  
WP5.SOP1 – Protocol for  
Brucella in vitro infection  
assay**

**Workpackage 5**

Responsible Partner: WBVR, ANSES

Contributing partners: APHA, BfR, FLI, IZSAM



## GENERAL INFORMATION

European Joint Programme full title	<i>Identification of <u>emerging Brucella</u> species: new threats for human and animals</i>
European Joint Programme acronym	One Health EJP JRP-17 IDEMBRU
Funding	This project has received funding from the European Union's Horizon 2020 research and innovation programme under Grant Agreement No 773830.
Grant Agreement	Grant agreement n° 773830
Start Date	01/01/2020
Duration	36 Months

## DOCUMENT MANAGEMENT

<b>Project SOP</b>	D-JRP17-WP5.SOP1 - Protocol for Brucella in vitro infection assay
<b>Project Acronym</b>	IDEMBRU
<b>Author</b>	Marielle van den Esker, Jeroen Koomen, Luca Freddi, Vitomir Djokic
<b>Other contributors</b>	Roland Ashford, Adrian Whatmore, Sascha Al Dahouk, Dirk Hofreuter, Daniela Prasse, Sandra Cavaco, Giuliano Garofolo, Flavio Sacchini, Fabrizio De Massis, Katuscia Zilli, Mihail Milanov, Hristo Daskalov, Ana Cristina Ferreira, Guillaume Girault, Acacia Ferreira Vicente, Claire Ponsart
<b>Due month of the report</b>	31 – Decembre 2022
<b>Actual submission month</b>	April 2022
<b>Type</b> <i>R: Document, report DEC: Websites, patent filings, videos, etc.; OTHER</i>	<i>OTHER</i> <b>Save date:</b> 26-Apr-23
<b>Dissemination level</b> <i>PU: Public (default) CO: confidential, only for members of the consortium (including the Commission Services)</i>	PU <b>This is the default setting.</b> If this project deliverable should be confidential, please add justification here (may be assessed by PMT): ..... .....
<b>Dissemination</b> <i>Author's suggestion to inform the following possible interested parties.</i>	OHEJP WP 1 <input type="checkbox"/> OHEJP WP 2 <input type="checkbox"/> OHEJP WP 3 <input type="checkbox"/> OHEJP WP 4 <input type="checkbox"/> OHEJP WP 5 <input type="checkbox"/> OHEJP WP 6 <input type="checkbox"/> OHEJP WP 7 <input type="checkbox"/> Project Management Team <input checked="" type="checkbox"/> Communication Team <input checked="" type="checkbox"/> Scientific Steering Board <input checked="" type="checkbox"/> National Stakeholders/Program Owners Committee <input checked="" type="checkbox"/> EFSA <input checked="" type="checkbox"/> ECDC <input checked="" type="checkbox"/> EEA <input checked="" type="checkbox"/> EMA <input checked="" type="checkbox"/> FAO <input checked="" type="checkbox"/> WHO <input checked="" type="checkbox"/> OIE <input checked="" type="checkbox"/> Other international stakeholder(s): ..... Social Media: ..... <b>Other recipient(s):</b> .....



# PROTOCOL FOR BRUCELLA IN VITRO INFECTION ASSAY

## Infection conditions

- MOI: 100
- Tested for cell line: HeLa, THP-1, BoMac, DH82
- Medium: RPMI + 10% FCS + 1%X NEAA/ DMEM + 10% FCS

### Day -4:

## Inoculation of bacteria on plates (BSL3)

- Streak *Brucella* on non-selective agar plates (eg TSA, HIS agar) from -80°C stocks (table below)
- Incubate plates at 37°C, 5% CO<sub>2</sub> and wait until colonies are visible

### Day -1

## Cell culture preparation: early morning

- 24 h before starting the experiment, trypsinize and homogenize cells from pre-culture
- Count cells (determine amount of dead cells)
- Seed cells in the following densities in order to obtain 95-100% confluence in each well:

Cell line	24-wells plate (cells/well)	Cells per cm <sup>2</sup>
HeLa	1 x 10 <sup>5</sup>	5.38 x 10 <sup>4</sup>
THP1	2 – 5 x 10 <sup>5</sup>	13.45 x 10 <sup>4</sup>

## Bacteria inoculation

Inoculate the *Brucella* strains in fresh liquid medium, so that 2\*10<sup>7</sup> bacteria per well are ready on the day of infection:

- Into a 25cm<sup>2</sup> flask add 15 mL non-selective medium (e.g. HIS broth or TSB)
- Check the purity of the *Brucella* plate and if ok, dissolve 1 to 2 colonies into the media
- Incubate the flask at +37°C, 5% CO<sub>2</sub>, shaking 170 rpm

Protocol is tested for the following reference strains:

NCTC number	Species	Biovar	Name	ATCC number
NCTC 10510	<i>Brucella suis</i>	biovar 2	Thomsen	
NCTC 10854	<i>Brucella canis</i>		RM-666	ATCC 23365
NCTC 10093	<i>Brucella abortus</i>	biovar 1	544	ATCC 23448
NCTC 10094	<i>Brucella melitensis</i>	biovar 1	16M	ATCC 23456
NCTC 10512	<i>Brucella ovis</i>		63/290	ATCC 25840
	<i>Brucella microti</i>		CCM 4915	
NCTC 10316	<i>Brucella suis</i>	biovar 1	1330	ATCC 23444;1330
NCTC 12890	<i>Brucella</i>	pinnipedialis		
NCTC 12891	<i>Brucella</i>	ceti		
	<i>Brucella canis</i>	Field strain		
BRU_2017_118	<i>Ochrobactrum intermedium</i>	X18004437-002	BSL3-HCU	Ringtestmonster Anses 2018
	<i>Brucella abortus</i>	Field strain		



## Day 0

### Infection of cell lines with *Brucella* (BSL3)

#### *Preparation of bacterial inoculum*

- In a 1.5 ml tube, put 1 ml of *Brucella* liquid culture (see day -1). If necessary, make more tubes for each strain
- Centrifuge at 4000 xg for 3 minutes
- Discard the supernatant and resuspend the pellet in 1 ml of PBS for washing
- Centrifuge at 4000 xg for 3 minutes
- Discard the supernatant and resuspend the pellet in 1 ml of PBS
- Measure OD600 and dilute inoculum with cell medium (RPMI + 10% FCS + 1X NEAA) to OD600 of 0.09 (OD/0.09 = dilution factor) in order to obtain a MOI of 100.
- Vortex inoculum
- Make serial dilutions (1:10) in RPMI + 10% FCS + 1%X NEAA) up to  $10^{-7}$ , plate last 4 dilutions to verify MOI
- 500µl of inoculum (*Brucella* in RPMI + 10% FCS + 1%X NEAA) is used per well for the infection

#### *Washing and infection of cells*

- Pipet medium off cells in plates that were seeded yesterday
- Wash the wells 2x with PBS
- Carefully remove PBS after last wash
- Add bacterial inoculum (500 ul) to each well
- Spin the plates 800 xg, 15 min --> timepoint 0 (start of infection)
- Incubate the 24 well plates for 1 hour at 37°C, 5% CO<sub>2</sub>
- Carefully remove bacterial inoculum
- Add 1 mL medium + 50 ug/ml gentamicin per well
- Incubate plates at 37°C, 5% CO<sub>2</sub> until sampling

## Timepoint of sampling

#### *Harvesting of cells for RNA isolation*

- After removal of the medium, immediately add 300 ul Trizol (see plate layout p 2)
- Incubate for 5 minutes at room temperature
- Scrape the cells loose with a pipet tip and put lysed cells in eppendorf tube
- Disinfect outside of tubes and take out of BSL3
- Isolate RNA according to protocol of kit supplier

#### *RT-PCR*

- Pipet off the medium of the wells
- Wash twice with PBS
- Pipet last PBS off very carefully
- Add 200 ul 0.1% - 0.2% TritonX-100 per well
- Incubate for 5 minutes at RT
- Homogenize cells with a P200 pipet tip
- Inactivate and lyse the cell suspension by pipetting 200 ul in 2 mL EasyMag lysis buffer
- Incubate overnight at RT
- Take out of BSL3 lab and proceed with DNA isolation followed by the PCR protocol provided

#### *CFU plating*

[Date]

[name of meeting] hosted by [name of host institution] in [location]



- Pipet medium out of well
- Wash the wells two times with PBS.
- Add 200  $\mu$ L of Triton X-100 (0.1/0.2% in PBS) for lysing the cells
- Scratch the cell with a p200 hooked.
- Homogenizing cell lysing and incubation 5 min.
- Serial dilution (in 96-wells plate, PBS) with 20  $\mu$ L in 180 $\mu$ L
- Plating 100 $\mu$ L of appropriate dilutions into TSA or BAB2 plates
- Incubation at 37°C +/- 5% CO<sub>2</sub> depending on Brucella species/strains

*Dilutions to be plated:*

<b>Time post infection (hours)</b>	<b>Brucella strain able to enters end replicates</b>	<b>Unknown Brucella strain</b>	<b>Brucella strain unable to enters end replicates</b>
1.5 – 2 (bacteria internalisation)	From 10 <sup>-2</sup> to 10 <sup>-4</sup>	From 10 <sup>0</sup> to 10 <sup>-4</sup>	From 10 <sup>0</sup> to 10 <sup>-2</sup>
5 – 6 (BCV maturation)	From 10 <sup>-2</sup> to 10 <sup>-4</sup>	From 10 <sup>0</sup> to 10 <sup>-4</sup>	From 10 <sup>0</sup> to 10 <sup>-2</sup>
20 – 24 (bacteria early replication)	From 10 <sup>-3</sup> to 10 <sup>-6</sup>	From 10 <sup>0</sup> to 10 <sup>-6</sup>	From 10 <sup>0</sup> to 10 <sup>-2</sup>
48 (end of bacteria replication)	From 10 <sup>-4</sup> to 10 <sup>-6</sup>	From 10 <sup>0</sup> to 10 <sup>-6</sup>	From 10 <sup>0</sup> to 10 <sup>-2</sup>